

VINDETO



ADL Dressing aid



www.vindetoprocurement.co.za

admin@vindeto.co.za 064 510 6179

DRESSING STICK



- Allows persons with limited mobility to dress without bending.
- Use the lightweight birch Dressing Stick with all types of clothing.
- The 16mm diameter dowel has a "C" hook at one end, and a push-pull hook at the other end.
- Plastic-coated hooks reduce scratches and snags.

DRESSING PAL™ SHOEHORN & DRESSING AID



- Multi-purpose dressing device combines a dressing stick and a shoehorn. Ideal for those who have limited reach or dexterity. Use push-pull end to put on clothing or as a reacher. Made of durable plastic with a vibrant purple foam handle. Available in both 24" and 30" lengths. Shoehorn tip measures 5.5" (14cm). Not made with natural rubber latex.

NORCO® BIG-GRIP™ BUTTON HOOK WITH ZIPPER PULL

- Get a secure, cushioned grip on hard-to-grasp buttons. The soft, cushioned grip makes the sometimes difficult task of buttoning clothes easier. Built with a Button Hook on one end of the handle, and a Zipper Pull on the opposite end. The built-up handle is made of a rubber-like material and features flexible ribbing that adapts to any grip. The firm 2" (5cm) long wire section is uniquely formed so it can slip around most buttons. Handle measures 4-1/4" (11cm) long and has a diameter of 1-1/2" (3.8cm). The total length measures 7.5" (19cm). Not made with natural rubber latex.



ZIP-GRIP ZIPPER PULL



- Snaps onto the zipper to extend the pull tab. The Zip-Grip Zipper Pull is ideal for persons who have difficulty grasping and pulling small zipper tabs. Just insert the peg of the Zip-Grip through the hole on a standard zipper, fold in half and press to snap shut. Ring diameter measures 1" (2.5cm). Sold in a package of 6

NORCO™ PLASTIC SHOE HORN

- Light weight classic shoehorn This durable, lightweight Norco™ Plastic Shoehorn can be grasped with ease. Measures 18" (45cm) long. Handle has a large hole for hanging.



NORCO™ PLASTIC SHOE HORN WITH HOOK

- Hooked handle pulls up socks. Norco™ Plastic Shoehorn with Hook offers a comfortable hold for fingers and may be used to pull up socks. Nylon loop for easy storage. Measures 16-1/4" (41cm) long.



EASY-SLIDE SHOE HORNS

- Easy-sliding shoehorns.
- Easy-Slide Shoehorns have an enameled steel shaft and rubber-like handle for easy grasping.
- They help guide feet into tight-fitting dress shoes, sport shoes or boots
- Available in four lengths.



NORCO® FOAM HANDLE SHOEHORNS



- Built-up foam handle is easy to grasp. The Foam Handle Shoehorns come with a looped string for a better grasp and positioning. Enameled steel shaft. Available in three lengths

BIG-GRIP™ SHOEHORNS



- Offers an extended reach and a secure, cushioned grip. Built-up handle made of a latex-free, rubber-like material with flexible ribbing that adapts to any grip. This handle is so unique, it carries a U.S. patent. Shaft is made of enameled steel. People who have limited mobility and difficulty bending over will appreciate the extra-long design. The textured surface prevents the foot from sticking to the shoehorns, helping it slide easily into the shoe. Handle diameter measures 1-1/2" (3.8cm).

FOOTFUNNEL™ SHOE ASSIST

- A great alternative to long handled shoe horns. The FootFunnel™ Shoe Assist is an ideal dressing aid for people who can't reach down. This device spreads the heel of the shoe, making it easy to guide the foot along the smooth surface and into the shoe



NORCO™ EASY-PULL™ SOCK AID



- Foot slides easily on the polyester lining while the sock remains secure. The flexible plastic core of the Norco™ Easy-Pull™ Sock Aid allows people to use just one hand to shape the trough easily to slide it into a sock. Large loop handles at the end of the strong, 31-1/4" (80cm) straps are great for persons with weak grasps. Polyester lining allows foot to slide easily along the trough. Terrycloth outer covering keeps the sock or stocking securely on the sock aid. Slits in the plastic core allow for contour bending around the heel. Compact and lightweight.

ACHIEVA™ SOCK-ASSIST™



- Unique features make this the most extraordinary sock aid on the market. The easy-to-use Achieva™ Sock-Assist™ lays flat on the floor, accommodates feet and calves of all sizes, has a non-stick surface and pulls the sock up the calf instead of releasing it at the ankle. Its innovative design allows easy donning of socks or knee-high nylons using one or two hands. Also works with elastic compression hose. Two side wings keep the sock on the Sock-Assist™ until it is released on the calf. Center tongue helps ease the sock onto the device. Flexible polypropylene “gives” to accommodate socks. Flat bottom surface prevents the Sock-Assist™ from rolling. Wide foot bed accommodates average to wide feet. Width ranges from 4.5" to 5.5". Length measures 7". Contoured side walls help hold the sock on the device. Rear flares prevent the cord from rubbing against the skin. Textured surface prevents skin from sticking to the sock aid. Available in two models, With Loop Handle and With Two Cord Handles. The 24" (61 cm) cord handles can be cut to adjust length.

NORCO® MOLDED SOCK AIDS



- A foam disc on the underside of the trough prevents the sock from sliding off before it is completely on the foot. Norco® Molded Sock Aids have a semi-flexible plastic trough that holds a sock or stockings in place while being pulled onto the foot. The trough may be adjusted using a heat gun. Good hand strength is required to flex the sock aid. Sock mold measures 9.6" x 5" x 3.2" (24 x 13 x 8cm). Available in two styles. The Norco® Molded Sock Aid with Single Cord Handle can be used with one or two hands. The 28" (71cm) cord can be trimmed to adjust the cord length. The Norco® Molded Sock Aid with Two Cord Handles features soft foam grips for comfortable use.

ENGLISH SOCK AID



- An economical solution for putting on socks. The English Sock Aid is lightweight and easy to use. The plastic trough can be bent as necessary to fit any size sock or nylon stocking. Indentations hold the sock securely until it is in place. The 32" (81cm) long cotton straps do not stretch during use. Straps can be adjusted as required by cutting or knotting. Tie the straps into a single loop if desired.

SOCK-EEZ®



- A foam disc on the underside of the trough prevents the sock from sliding off before it is completely on the foot. Norco® Molded Sock Aids have a semi-flexible plastic trough that holds a sock or stockings in place while being pulled onto the foot. The trough may be adjusted using a heat gun. Good hand strength is required to flex the sock aid. Sock mold measures 9.6" x 5" x 3.2" (24 x 13 x 8cm). Available in two styles. The Norco® Molded Sock Aid with Single Cord Handle can be used with one or two hands. The 28" (71cm) cord can be trimmed to adjust the cord length. The Norco® Molded Sock Aid with Two Cord Handles features soft foam grips for comfortable use

HEEL GUIDE™ COMPRESSION STOCKING AID



- Eases donning of compression hosiery. The unique indented channel on the center cone acts as a heel guide and keeps the foot in position while pulling on the compression hose. Heel Guide™ Compression Stocking Aid has an extra wide cone is ideal for people with edema or larger feet. Sturdy 17" (43cm) high handles reduce the need to bend or stretch

COILERS™ SHOELACES



- Easy-to-use, adjustable Coilers™ never need tying. Coilers™ Shoelaces may be tightened or loosened, but never need tying. They have no pressure points and allow natural movement. Fits up to eight pairs of eyelets. Sold in pairs.

NORCO™ ELASTIC SHOELACES



- Turn a lace-up shoe into a slip-on shoe. Norco™ Elastic Shoelaces are available in two varieties, Regular and Deluxe. Lace and tie shoes in the normal manner just once. Shoes can be slipped on or off easily without tying and untying. Regular measures 1/8" (3.2mm) wide. Deluxe measures 3/8" (9.5mm) wide. Available in two lengths: 24" (61cm) or 37" (94cm), and three colors: Black, Brown or White

EZ REACHER®



Locking EZ Reacher®



- Durable reacher design stands up to everyday use.
- The EZ Reacher® relieves the strain of bending, stretching or stooping to pick up items beyond reach.
- Lightweight aluminum jaws with flexible, non-slip gripping cups conform easily to the shape of any item.
- Jaws open to 8.9cm wide.
- Made of lightweight rustproof aluminum with an easy-to-squeeze trigger.
- Handle color may vary.
- Lifetime warranty.
- Available in three models: Short, Standard, and Long.
- The Locking EZ Reacher® has all the features listed above, but maintains a hold on items without constant trigger pressure.
- Pull down the lever to keep the jaws locked.
- Release the lock by fully retracting the trigger and lifting the lever up.