

VINDETO



Mobility Aid



www.vindetoprocurement.co.za

admin@vindeto.co.za 064 510 6179

ADJUSTABLE RECIPROCAL FOLDING TWO STEP WALKER



Dimensions

Length	53cm
Width	47cm
Height	73 - 85cm (Adjustable)

Weight Limit

Weight Limit	110kg
--------------	-------

- An easy-to-use reciprocal folding walker with stand assist. The reciprocal walking frame walks with the user. There is no need to lift the walker forward with each step, rather this walker when used correctly, steps with the user. This walker may also be used as a rigid frame.
- This walker can be folded with the push of a quick single release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.



FOLDING WALKER WITH WHEELS

- An easy-to-use folding walker designed as an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick double release button allowing for efficient storage and travel. The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.
- Wheels can be added to convert the walker into a rollator designed to be used primarily by people who are recuperating from leg or back injuries as well as hip and knee replacements. The two wheels at the front make for easier and more energy efficient propulsion whilst still providing a stable base of support.
- Two push button release system
- Also available as a rollator – wheels sold separately
- Weight limit Adult 110kg, Youth 80kg, Paediatric 60kg
- Height from handle to ground Adult 80-87cm, Youth 66-75cm, Paediatric 56-63cm



FOLDING WALKER – ONE PUSH BUTTON



Dimensions

Length	58cm
Width	49cm
Height	79 - 97cm (Adjustable)

Weight Limit

Weight Limit	110kg
--------------	-------



- An easy-to-use folding walker designed as an ideal and dependable solution for those who require a supportive yet durable aid for their daily activities. This walker can be folded with the push of a quick single release button allowing for efficient storage and travel.
- The lightweight aluminium frame features soft grip handles and height adjustable legs to ensure that the set-up is best suited to the user's anatomy.
- Available in natural colour only
- Length 58cm, width 49cm, height 79-97cm
- Weight limit 110kg

FOLDING WALKER TWO PUSH BUTTON

- Designed to be used primarily by people who are unable to take weight through their wrists and hands, for example, people with Rheumatoid Arthritis of the hands who cannot grasp handles of frames.
- Two push button release system
- Also available as a rollator – wheels sold separately
- The gutters are detachable
- Weight limit 110kg



Dimensions

Length	58cm
Width	49cm
Height	79 - 97cm (Adjustable)

Weight Limit

Weight Limit	110kg
--------------	-------

FOLDING WALKER (YOUTH)



FOLDING WALKER STANDARD KIDDIES EXTRA SMALL

- Length 47cm, width 58cm, height 58-68cm

- Length 35cm, width 39cm, height 49-54cm

ALUMINUM LIGHT WEIGHT ROLLATOR WITH OPEN BASKET



Features:

- Ergonomic Handgrips
- Light touch hand brake.
- With a sling backrest and mesh seat.
- Height adjustable handles
- With a strong aluminum basket for storage
- Large wheels of 15 cm
- Padded seat 35 cm x 30cm
- Can be folded into a compact size for packing
- Curved backrest for a comfortable place to sit



2 models available- 6 inches wheels and 8 inches wheels.

STANDARD ELBOW CRUTCH (ADULT)



- An ergonomically designed elbow crutch for adult patients and primarily to assist with recovery, improve mobility and provide support for lower extremity injuries.
- Available colours are natural, black, blue and red
- Sold in pairs
- Extensions Available (55cm)
- Weight limit Adult 120kg, Youth 80kg and Paediatric 60kg
- Height range from handle to ground Adult 70,3-92,7cm, Youth 51-73cm and Paediatric 41-63cm

PLASTIC TOP ELBOW CRUTCH

- An ergonomically designed elbow crutch for adult patients and primarily to assist with recovery, improve mobility and provide support for lower extremity injuries.
- Sold in pairs
- Glass fitted polypropylene top (forearm support) for extra strength
- Contoured arm cuffs and grips moulded for comfort and better stability
- Single pin for adjustable height
- Extensions available (55cm)



- Available colours for Adult black, dark blue, aqua, grey, green, orange, yellow, red, and burgundy,
- Weight limit Adult 120kg,
- Height range from handle to ground Adult 70-90cm,



- Youth and Paediatric aqua, red and yellow
- Weight limit Youth 80kg, Paediatric 60kg
- Height range from handle to ground Youth 52-74,5cm, Paediatric 42-64cm

SOLID PERMANENT CRUTCH

- A permanent crutch Ideal for patients who require long-term walking support.
- Available colours are natural, black, blue and red
- Sold in pairs
- Length can be cut to size
- Standard option weight limit 120kg
- Heavy duty option weight limit 150kg
- Height range from cuff to ground 120cm






ADJUSTABLE CLASSIC WALKING STICK

- The perfect solution for the elderly, the physically challenged, or patients who have trouble walking or standing. The design is suitable for left and right handed users.
- Available colours are white, natural, black and bronze
- Weight limit 110kg
- Adjustable height from handle to ground 70 – 90cm

CANE WALKING STICK

- A cane walking stick which is primarily used to aid walking, provide postural stability and support, or assist in maintaining a good posture.
- Available in a natural colour with standard colours available on request
- Weight limit 100kg
- Height range from handle to ground 72–94cm




MWS-CH-1351




MWS-??-1352



MSB-A-1354



MC-??-1353



MWC-CH-1354

TRIPODS & QUADRUPED

A tetrapod is another name for a quad cane – a walking stick with four legs. Tripod and tetrapod walking sticks provide sturdy, lightweight alternatives to conventional sticks and canes. These items significantly improve stability for those with mobility or balance issues. As well as offering improved support, ...



WHAT IS THE DIFFERENCE BETWEEN A WALKER & ROLLATOR?

Walkers are either without wheels or have wheels on the front two legs, while rollators have wheels on all four legs (three or four, depending on the model). The most supportive option is a regular walker without wheels, but you must lift it with each stride. Rollators offer support without slowing down your walking pace





GUTTER TOP

- Sold in pairs
- Soft sponge on base plate and Velcro strap
- Adjustable plastic handle
- Available in natural colour only
- Weight limit 110kg



ROLLATOR WHEEL SET

Used with the folding walker to create a rollator.

- Available in natural anodised colour only
- Pack Size: Set / 2 Units
- Wheel Size: 5" / 125 mm



PLASTIC CUFFS

- Available in adult, youth and paediatric sizes in black only.



STANDARD FOOT

Standard crutch foot made from black PVC.

- Available in 16, 22 and 24mm diameter
- 10 units per pack
- Easy to attach
- Should be replaced regularly
- Thick base and galvanised washer inside foot to help prevent splitting



SUPER SOFT GRIPS

Available in adult size and made from thick black rubber.